

CARING FOR SENSORY DEVICES

- Check battery with battery tester
- Check earmold and tubing for cracks & clogs and moisture
- Check case for damage
- Check switches for mobility
- Check volume setting
- Check for feedback (“eee”)
 - If present, check connections, check for cracks in tubing, check fit of ear mold
- 6 Sound Test (ah, oo, ee, sh, s, m)

CARING FOR HEARING AIDS

- Keep away from heat and moisture
- Store batteries and aids away from small children
- Replace dead batteries immediately
- Clean hearing aids as instructed
- Keep hairspray, gel, etc. away from aids
- Turn off hearing aids when not in use

CARING FOR COCHLEAR IMPLANTS

- Keep dry
- Avoid activities that could cause severe blow to the head; use protective head gear
- Avoid exposure to static electricity (plastic slides, plastic ball rooms, unprotected computer monitors)

Test cell phones before purchase for interference